

Product datasheet for AP08095PU-S

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Phospholipase C gamma 1 (PLCG1) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: IHC, WB

Recommended Dilution: Western Blot (1/500-1/1000).

Immunohistochemistry (1/50-1/100).

Reactivity: Human, Mouse, Rat

Host: Rabbit

Clonality: Polyclonal

Immunogen: The antiserum was produced against synthesized non-phosphopeptide derived from human

PLC-gamma 1 around the phosphorylation site of Tyrosine 771 (P-D-YP-G-A).

Specificity: This antibody detects endogenous levels of total PLC-gamma 1 protein.

Formulation: PBS (without Mg2+ and Ca2+), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol.

State: Aff - Purified

State: Liquid purified Ig fraction.

Concentration: lot specific

Purification: Immunoaffinity Chromatography.

Conjugation: Unconjugated

Storage: Store the antibody (in aliquots) at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: One year from despatch.

Gene Name: phospholipase C gamma 1

Database Link: Entrez Gene 5335 Human

P19174

Synonyms: PLC1, Phospholipase C-gamma-1, PLC-gamma-1, PLC gamma1 Phospholipase C-II, PLC-II, PLC-

148





Product images:




Figure 2. Western blot analysis of extracts from 293 cell using PLC-gamma 1 Antibody AP08095PU (Lane 1 and 2).




Figure 1. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using PLC-gamma 1 Antibody AP08095PU.