

Product datasheet for **AP08138FC-N**

Turkey IgG (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Turkey IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Immunofluorescence: < / = 1 µg/10e6 cells. ELISA. Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry.
Reactivity:	Turkey
Host:	Goat
Isotype:	IgG
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified IgG fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography on Chicken IgG covalently linked to agarose.
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)