

Product datasheet for **AP08138HR-N**

Turkey IgG (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Turkey IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/4,000-1/8,000.
Reactivity:	Turkey
Host:	Goat
Isotype:	IgG
Formulation:	50% PBS/50% Glycerol, pH 7.4 without preservatives. Label: HRP State: Liquid purified IgG fraction. Label: Horseradish Peroxidase
Purification:	Affinity Chromatography on Turkey IgY covalently linked to agarose
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Source: Pooled antisera from goats hyperimmunized with IgY.



[View online »](#)