

Product datasheet for **AP08138PU-N**

Turkey IgG (H+L chain) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Turkey IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Immunofluorescence. ELISA. Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry.
Reactivity:	Turkey
Host:	Goat
Isotype:	IgG
Formulation:	100 mM Borate Buffered Ssaline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified IgG fraction.
Concentration:	lot specific
Purification:	Affinity Chromatography on Chicken IgG covalently linked to agarose.
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)