

## Product datasheet for **AP08138RP-N**

### Turkey IgG (H+L chain) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Turkey IgG (H+L chain) Goat Polyclonal Antibody
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	<b>Immunofluorescence.</b> <b>ELISA.</b> <b>Western blotting.</b> <b>Dot- and Slot-Immunoblotting.</b> <b>Immunohistochemistry.</b>
Reactivity:	Turkey
Host:	Goat
Isotype:	IgG
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified IgG fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Purification:	Affinity Chromatography on Chicken IgG covalently linked to agarose.
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)