

Product datasheet for **AP09111AG-N**

Protein A Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Protein A [Staphylococcus aureus]
Specificity:	This antibody reacts to Protein A [Staphylococcus aureus].
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide Label: Agarose State: Liquid IgG fraction
Conjugation:	Agarose
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.



[View online »](#)