

## Product datasheet for AP09121TR-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Protein A Goat Polyclonal Antibody**

**Product data:** 

Product Type: Primary Antibodies

Applications: ELISA, FC, IF, IHC

**Recommended Dilution:** Suitable for immunomicroscopy and flow cytometry or FACS analysis as well as other

antibody based fluorescent assays requiring lot-to-lot consistency.

Flow cytometry: 1:2,000 - 1:10,000. Immunoflourescence: 1:500 - 1:2,500.

Host: Goat lgG

**Clonality:** Polyclonal

Immunogen: Protein A [Staphylococcus aureus]
Specificity: This antibody detects Protein A.

Assay by immunoelectrophoresis resulted in a single precipitin arc against purified and partially purified Protein A [Staphylococcus aureus], anti-Goat IgG and anti-Goat Serum.

Formulation: 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 as buffer, 10 mg/ml Bovine

Serum Albumin (BSA) IgG and Protease free as stabilizer, 0.01% (w/v) Sodium Azide as

preservative Label: Texas Red

State: Lyophilized Ig fraction

Label: (TM) Sulfonyl Chloride (Molecular Weight 625 daltons)

Absorption emission: 596 nm / 620 nm

Molar radio: 1.5 moles Texas Red(TM) per mole of Goat IgG

**Reconstitution Method:** Restore with 2.0 ml of deionized water (or equivalent).

**Concentration:** lot specific

**Purification:** Multi-step process which includes delipidation, salt fractionation and ion exchange

chromatography followed by extensive dialysis against the buffer

Conjugation: Texas Red





## Protein A Goat Polyclonal Antibody – AP09121TR-N

Storage: Store vial at 2-8° C prior to restoration. For extended storage aliquot contents and freeze at -

20° C or below. Avoid cycles of freezing and thawing. Centrifuge product if not completely clear after standing at room temperature. This product is stable for two weeks at 2-8° C as an

undiluted liquid. Dilute only prior to immediate use.

**Stability:** Shelf life: One year from despatch.