

Product datasheet for **AP09205TC-N**

GST-Tag Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	FC
Recommended Dilution:	Immunofluorescence. Flow Cytometry.
Reactivity:	Schistosoma japonicum
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	The immunogen is full length GST isolated from Schistosoma japonicum.
Specificity:	This antibody reacts to Glutathione-S-Transferase (GST).
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) as stabilizer and 0.01% (w/v) Sodium Azide as preservative Label: TRITC State: Lyophilized Label: Tetramethylrhodamine isothiocyanate Absorption emission: 550 nm / 570 nm Molar ratio: 2.8 moles TRITC per mole of Goat IgG
Reconstitution Method:	Restore with 0.1 ml of deionized water (or equivalent).
Concentration:	lot specific
Conjugation:	TRITC
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P08515



[View online »](#)

- Background:** GST (Glutathione S-Transferase) is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.
- Synonyms:** Glutathione S-transferase Tag, GST26-Tag