

Product datasheet for AP09230HR-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

HA Epitope Tag (YPYDVPDYA) Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA, IHC, WB Recommended Dilution: ELISA: 1/15,000.

Western Blot: 1/1,000-1/5,000.

Immunohistochemistry: 1/500-1/2,000.

Note: Anti-HA is optimally suited for monitoring the expression of HA-tagged fusion proteins. As such, anti-HA/HA can be used to identify fusion proteins containing the HA epitope. The antibody recognizes the HA epitope tag fused to the amino- or carboxy- termini

of targeted proteins, as expressed in many commonly used expression vectors. This

antibody has been tested by ELISA and western blotting against both the immunizing peptide

and HA containing recombinant proteins. Although not tested, this antibody is likely

functional for other immunodetection techniques.

Host: Rabbit Isotype: IgG

Clonality: Polyclonal

Immunogen: Purified from whole Rabbit serum prepared by repeated immunizations with the 9-aa epitope

tag peptide YPYDVPDYA (114-122) from Hemagglutinin Influenza conjugated to KLH using maleimide. A residue of cysteine was added to the carboxy terminal end to facilitate coupling

Specificity: This antibody is directed against the HA motif and is useful in determining its presence in

various assays. This polyclonal anti-HA tag antibody detects over-expressed proteins

containing the HA epitope tag.

In Western blotting of bacterial extracts, the antibody does not cross-react with endogenous

oroteins.

Affinity purification of the polyclonal antibody results in very low background levels in assays

and low cross-reactivity with other cellular proteins.

Formulation: 0.02M Potassium Phosphate, 0.15M Sodium Chloride, pH 7.2

Label: HRP

State: Lyophilized purified IgG fraction

Stabilizer: 10 mg/ml BSA

Preservative: 0.01% (w/v) Gentamicin Sulfate

Label: Horseradish Peroxidase





HA Epitope Tag (YPYDVPDYA) Rabbit Polyclonal Antibody - AP09230HR-N

Reconstitution Method: Restore with 0.1 ml of deionized water (or equivalent).

Concentration: lot specific

Purification: Affinity Chromatography

Conjugation: HRP

Storage: Prior to reconstitution store at 2-8°C.

Following reconstitution store the antibody undiluted at 2-8°C for one month

or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Background: Epitope tags are short peptide sequences that are easily recognized by tag-specific

antibodies. Due to their small size, epitope tags do not affect the biochemical properties of the tagged protein. Most often, sequences encoding the epitope tag are included with the target DNA at the time of cloning to produce fusion proteins containing the epitope tag sequence. This allows anti-epitope tag antibodies to serve as universal detection reagents for any tag containing protein produced by recombinant means. This means that anti-epitope

tag antibodies are a useful alternative to generating specific antibodies to identify,

immunoprecipitate or immunoaffinity purify a recombinant protein. The anti-epitope tag antibody is usually functional in a variety of antibody-dependent experimental procedures. Expression vectors producing epitope tag fusion proteins are available for a variety of host

expression systems including bacteria, yeast, insect and mammalian cells.