

## Product datasheet for **AP09240PU-N**

### **Btb1 (1205-1217) Rabbit Polyclonal Antibody**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Applications:</b>	ELISA
<b>Recommended Dilution:</b>	ELISA: 1:15,000 to 1:65,000.
<b>Reactivity:</b>	Human
<b>Host:</b>	Rabbit
<b>Isotype:</b>	IgG
<b>Clonality:</b>	Polyclonal
<b>Immunogen:</b>	Synthetic peptide corresponding to aa 1205-1217 of Human Btb1
<b>Specificity:</b>	This antibody detects Btb1.
<b>Formulation:</b>	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.01% (w/v) Sodium Azide State: Aff - Purified State: Liquid sterile filtered Ig fraction
<b>Concentration:</b>	lot specific
<b>Purification:</b>	Affinity chromatography
<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody at 2 - 8 °C up to one month or (in aliquots) at -20 °C for longer. Avoid repeated freezing and thawing. Should this product contain a precipitate we recommend microcentrifugation before use.
<b>Stability:</b>	Shelf life: one year from despatch.
<b>Background:</b>	Btb1 is a currently uncharacterized human ankyrin repeat protein BAA92655. This protein has been suggested to bind to CUL3.



[View online »](#)