

Product datasheet for **AP09317PU-N**

REC107 (26-35) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/5,000 - 1/25,000. Western Blot: 1/1,000 - 1/10,000.
Reactivity:	Saccharomyces cerevisiae
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Synthetic peptide corresponding to amino acids 26-35 of Saccharomyces cerevisiae Mer2 protein
Specificity:	This antibody is directed against the Saccharomyces cerevisiae Mer2 protein.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 0.01% (w/v) Sodium Azide State: Aff - Purified State: Liquid purified Ig
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P21651




[View online »](#)

Background:

Mer2 (also known as meiotic recombination 2 protein) is a chromosomal protein that is critical for meiotic recombination and progression. It is phosphorylated at two serine residues, S30 and S271 by the yeast Cdk1 cyclin- dependent kinase homolog. This phosphorylation is S-phase specific, and thus has the potential to be a specific assay for S-phase cyclin-dependent kinases. Moreover, there are hints that the phosphorylation may be a mark of replication fork passage, which would indicate that Sphase CDK associates with the replication fork.

Synonyms:

Meiotic recombination 2 protein, REC107, YJR021C, J1462

Product images:

Western blot using affinity purified anti-*S.cerevisiae* Mer2 antibody shows detection of phosphorylated and unphosphorylated Mer2 in wild type, phosphatase treated and mutant cells. Lane 1 contains Mer2-myc protein detected in wild type cells after first immunoprecipitating the protein using anti-myc antibody. Cells were harvested 4 h after the initiation of meiosis and therefore contain mostly phosphorylated Mer2. Lane 2 contains the same preparation after treatment with phosphatase. Lane 3 contains Mer2-S30A protein as a phosphorylation control. This antibody is reactive with both phosphorylated and unphosphorylated Mer2 at the S30 position. The primary antibody was used at a 1:5,000 dilution.