

Product datasheet for **AP09429HR-N**

Biotin (F(ab')₂ fragment) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Western Blot: 1/200 - 1/1,000. Dot Blot. ELISA: 1/2,000 - 1/10,000. Immunohistochemistry: 1/100 - 1/500.
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Biotin conjugated to Keyhole Limpet Hemocyanin (KLH)
Specificity:	This is a Peroxidase conjugated F(ab) ₂ fragment of Anti-Biotin.
Formulation:	0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2 containing 10 mg/ml Bovine Serum Albumin (BSA) and 0.01% (w/v) Sodium Azide Label: HRP State: Lyophilized F(ab) ₂ fragment Label: Horseradish Peroxidase (Molecular Weight 40,000 daltons)
Reconstitution Method:	Restore with 1.0 ml of deionized water (or equivalent). For extended storage reconstitute product with 50% glycerol instead of water
Concentration:	lot specific
Purification:	Immunoaffinity chromatography
Conjugation:	HRP
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)

Background: Biotin is a water soluble vitamin, generally classified as a B complex vitamin, also called vitamin B4. After the initial discovery of biotin, nearly forty years of research were required to establish it as a vitamin. Biotin is required by all organisms but can only be synthesized by bacteria, yeasts, molds, algae, and some plant species. Biotin is required as prosthetic group of enzymes involved in incorporation of carbon dioxide into organic compounds. Biotin has a MW of 244 Da.

Synonyms: Vitamin B7, Vitamin H