

Product datasheet for AP09487PU-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cyclin B1 (CCNB1) pSer147 Rabbit Polyclonal Antibody

Product data:

Clonality:

Product Type: Primary Antibodies

Applications: IF, IHC, WB

Recommended Dilution: Immunofluorescence: 1/100 - 1/200.

Polyclonal

Immunohistochemistry on Paraffin Sections: 1/50 - 1/100.

Western Blot: 1/500 - 1/1000; Incubate membrane with diluted antibody in 5% nonfat milk,

1X TBS, 0,1% Tween-20 at 4°C with gentle shaking, overnight.

Reactivity: Human

Host: Rabbit

Immunogen: Synthetic phosphopeptide derived from human Cyclin B1 around the phosphorylation site of

Serine 147 (A-F-Sp-D-V).

Specificity: This Antibody detects endogenous levels of Cyclin B1 only when phosphorylated at Serine

147.

Formulation: PBS (without Mg2+ and Ca2+), pH 7.4, 150 mM NaCl, 0.02% Sodium Azide and 50% Glycerol.

State: Aff - Purified

State: Liquid purified Ig fraction

Concentration: lot specific

Purification: Immunoaffinity Chromatography

Conjugation: Unconjugated

Storage: Store the antibody (in aliquots) at -20°C.

Avoid repeated freezing and thawing.

Stability: Shelf life: One year from despatch.

Gene Name: cyclin B1

Database Link: Entrez Gene 891 Human

P14635





Background:


In mammals, cyclin B associates with inactive p34 cdc2 which facilitates phosphorylation of p34 cdc2 at amino acids 14 Thr and 15 Tyr. This maintains the inactive state until the end of G2 phase. The inactive cyclin B p34 cdc2 complex continues to accumulate in the cytoplasm until the completion of DNA synthesis, when Cdc25, a specific protein phosphatase, dephosphorylates amino acids 14Thr and 15Tyr of p34 cdc2, rendering the complex active at the G2 / M boundary. This mitotic kinase complex remains active until the metaphase / anaphase transition when cyclin B is degraded. This degradation process is ubiquitin dependent and is necessary for the cell to exit mitosis. Therefore, cyclin B p34 cdc2 plays a critical role in G2 to M transition. Two alternative transcripts have been found, a constitutively expressed transcript, and a cell cycle-regulated transcript that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites.

Synonyms: CCNB1, CCNB, Cyclin-B1

Protein Families: Druggable Genome, Stem cell - Pluripotency

Protein Pathways: Cell cycle, Oocyte meiosis, p53 signaling pathway, Progesterone-mediated oocyte maturation

Product images:



Western Blot analysis of extracts from MDA cells untreated or treated with Anisomycin using Cyclin B1 (pSer147) antibody






Figure 1. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Cyclin B1 pSer147 Antibody.




Figure 2. Immunofluorescence staining of methanol-fixed HeLa cells using Cyclin B1 pSer147 Antibody (Red).