

Product datasheet for **AP09542SU-N**

Periplakin (PPL) Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, IP, WB
Recommended Dilution:	Immunoprecipitation. Immunofluorescence. Immunohistochemistry on Frozen Sections: 1/50- 1/100. Western Blot: 1/1000-1/2000 (ECL). <i>Incubation Time:</i> 1 h at RT.
Reactivity:	Bovine, Human
Host:	Guinea Pig
Clonality:	Polyclonal
Immunogen:	Synthetic peptides of human periplakin (P1 a/b, P2a/b, P3a/b), coupled to KLH
Specificity:	Specific for Periplakin, a 190 kD protein localized e.g. in epidermis (cornified envelope) and cortical lens fiber cells.
Formulation:	State: Serum State: Liquid Stabilized Antiserum Preservative: 0.09% Sodium Azide
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Stability:	Shelf life: one year from despatch.
Gene Name:	periplakin
Database Link:	Entrez Gene 5493 Human O60437
Synonyms:	KIAA0568



[View online »](#)