

## Product datasheet for **AP09638PU-N**

### 17-alpha-Boldenone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 0.25 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	17-alpha-Boldenone(3)-BTG
Specificity:	This antibody recognizes 17-alpha-Boldenone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Boldenone is an anabolic steroid developed for veterinary use, mostly for treatment of horses. As it is not designed to be used by humans, it is only available at veterinary clinics.



[View online »](#)