

Product datasheet for **AP09643PU-N**

17-alpha-Trenbolone Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	17alpha-Trenbolone (3)-BTG
Specificity:	This antibody reacts to 17alpha-Trenbolone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Trenbolone is an anabolic steroid used mainly as a growth substance in animals. It is one of the best effective anabolic compounds, promoting protein synthesis, as well as creating a positive nitrogen balance.



[View online »](#)