

## Product datasheet for **AP09644PU-N**

### 17-beta-Boldenone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	17beta - Boldenone (3)-BTG
Specificity:	This antibody reacts to 17beta-Boldenone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Boldenone is an anabolic steroid developed for veterinary use, mostly for treatment of horses. As it is not designed to be used by humans, it is only available at veterinary clinics.



[View online »](#)