

## Product datasheet for **AP09645PU-L**

### 17-beta-Estradiol (-Long Linker) Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Estradiol (6)-short linker-BSA
Specificity:	This antibody reacts to 17beta-Estradiol (-Long Linker).
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	17 beta Estradiol is a potent mammalian estrogenic hormone that is produced in the ovary, placenta, testis, and possibly the adrenal cortex.



[View online »](#)