

## Product datasheet for **AP09649PU-N**

### 17-beta-19-Nor-Testosterone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 10µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	17beta -19- Nortestosterone(3)-BTG
Specificity:	This antibody reacts to 17beta-19-Nortestosterone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	19 Nortestosterone is a steroid with androgenic as well as anabolic properties that occurs normally in the human body. It is also known as the popular anabolic steroid Nandrolone.



[View online »](#)