

Product datasheet for **AP09659PU-L**

Albendazole Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Albendazole (5-propylthio)-BTG
Specificity:	This antibody reacts to Albendazole.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Albendazole, methyl N-(6-propylsulfanyl-1H-benzimidazol-2-yl)carbamate, is a benzimidazole broad-spectrum anthelmintic effective against many diseases caused by parasitic worms. It is used therapeutically in the treatment of Helminthiasis in man and animal.
Synonyms:	Albenza, Bilutac, Eskazole, Proftril, Valbazen, Zental, Zentel



[View online »](#)