

## Product datasheet for **AP09661PU-L**

### 4-Nitrophenyl AHD Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 5.0µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Aminohydantoin-BTG
Specificity:	This antibody reacts to 4-Nitrophenyl AHD.
Formulation:	20mM Phosphate, 150mM Sodium Chloride pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Nitrofurantoin, Furadantin, Macrochantin, Macrobid, Furalan, Nitrofuradantin, Macrochantina, Chemiofuran, Cistofuran, Furadantine



[View online »](#)