

Product datasheet for **AP09672PU-L**

Buprenorphine Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 11.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Buprenorphine (3)-BSA
Specificity:	This antibody reacts to Buprenorphine. Cross-reactivity: Buprenorphine 100% Morphine 1.2% Morphine-3-glucuronide 1.2% 6-Monoacetylmorphine 1.2% Codeine 1.2% Oxycodone 1.2% Noroxycodone 1.2% Hydromorphone 1.2% Hydrocodone 1.2% Duhydrocodone 1.2% Thebaine 1.2%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)