

## Product datasheet for **AP09675PU-L**

### 4-Nitrophenyl AMOZ Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 0.625-10 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Carboxybenzyl-AMOZ-BTG
Specificity:	Recognizes 4-NITROPHENYL-AMOZ
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Nitrofurantoin, Furadantin, Macrochantin, Macrobid, Furalan, Nitrofuradantin, Macrochantina, Chemiofuran, Cistofuran, Furadantine



[View online »](#)