

Product datasheet for **AP09680PU-N**

Chlormadinone Acetate Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Chlormadinone acetate (3)-BTG
Specificity:	This antibody recognizes Chlormadinone Acetate.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Chlormadinone acetate is a progestin used in oral contraceptive products.
Synonyms:	Gestafortin, Matrol, Menstridyl, 17alpha-Acetoxy-6-chloro-4,6-pregnadiene-3,20-dione, 6-Chloro-6-dehydro-17alpha-acetoxyprogesterone



[View online »](#)