

Product datasheet for **AP09687PU-L**

Codeine Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Codeine (6)-BSA
Specificity:	This antibody reacts to Codeine.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Codeine is an alkaloid narcotic, C ₁₈ H ₂₁ NO ₃ , derived from opium or morphine and used as a cough suppressant, analgesic, and hypnotic.
Synonyms:	Methylmorphine, Coducept, l-Codeine, Isocodeine, O3-Methylmorphine, Codeina, Codeinum



[View online »](#)