

## Product datasheet for **AP09700PU-L**

### Dehydroepiandrosterone / DHEA Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.625 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	DHEA-BTG
Specificity:	This antibody recognizes Dehydroepiandrosterone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Dehydroepiandrosterone (DHEA) is a major secretory product of the human adrenal cortex during intrauterine development as well as during adulthood. DHEA circulates in the bloodstream as DHEA-sulfate (DHEAS). DHEA has been implicated in a wide range of physiological roles including aging, immunology, learning and memory and obesity.
Synonyms:	3-beta Hydroxyandrost 5 en 17 one, 3-beta hydroxy 5 androsten 17 one, Dehydroisoandrosterone, Prasterone, Trans Dehydroandrosterone



[View online »](#)