

Product datasheet for **AP09703PU-L**

Dienestrol (p-OH)-BTG Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.078µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Dienestrol (p-OH)-BTG
Specificity:	This antibody reacts to Dienestrol.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Dienestrol is a synthetic, non-steroidal estrogen. It is commonly sold as a cream for intra-vaginal use.
Synonyms:	3, 4-Bis(4-hydroxyphenyl)-2, 4-hexadiene



[View online »](#)