

## Product datasheet for **AP09712PU-N**

### Epiandrosterone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1.25 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Epiandrosterone (17)-BTG
Specificity:	This antibody recognizes Epiandrosterone.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	3-epiandrosterone, 3-beta-androsterone, 3-beta-hydroxyetioallocholan-17-one, 5-alpha-androstan-17-one, 3-beta-hydroxy-, 5?-androstan-17-one, 3-beta-hydroxyetiollocholan-17-one



[View online »](#)