

Product datasheet for **AP09715PU-L**

Estriol Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.325 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Estriol (6)-BTG
Specificity:	This antibody recognizes Estriol (6).
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.4 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Estriol is a hydroxylated metabolite of ESTRADIOL or ESTRONE that has a hydroxyl group at C3-beta, 16-alpha, and 17-beta position. Estriol is the main estrogen in the blood and urine of pregnant females. During pregnancy, large amount of estriol is produced by the placenta. Serum estriol determinations reflect the status of the fetus and are used to evaluate fetal distress and placental function. Estriol increases until term, reaches a plateau at the 40th week of gestation and then declines.
Synonyms:	Oestriol
Note:	Vial should be centrifuged briefly before opening to ensure all materials is removed from the cap.



[View online »](#)