

Product datasheet for **AP09718PU-N**

Ethinylestradiol Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Ethinylestradiol(3)-BTG
Specificity:	This antibody recognizes Ethinylestradiol.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Ethinylestradiol is a semisynthetic alkylated Estradiol with a 17-alpha-ethinyl substitution. It has high estrogenic potency when administered orally, and is often used as the estrogenic component in oral contraceptives.
Synonyms:	Estinyl, Ethinyl Estradiol, Ethinylestradiol, Etinestrol, Etinestryl, Etinoestryl, Etistradiol, Follicoral, Ginestren
Note:	Chemical Formula: C ₂₀ H ₂₄ O ₂ Mol. Mass: 296.40336 g/mol



[View online »](#)