

Product datasheet for **AP09725PU-N**

Flunixin (COOH) Sheep Polyclonal Antibody

Product data:


Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.625 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Flunixin(COOH)-BTG
Specificity:	This antibody recognizes Flunixin
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative. State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Flunixin meglumine is a member of the nonsteroidal anti-inflammatory drugs (NSAID's) used for pain relief (particularly colic), counteract the effects of trauma, bring down fever, and block the effects of endotoxins on tissue. Frequently used to treat medical colics though care must be exercised as it may prevent early recognition of surgical problems by covering up symptoms. Flunixin is useful as a medium duration analgesic and anti-inflammatory drug. Like other NSAIDs, it also has some anti-pyretic properties. It is used for visceral or musculoskeletal pain in horses.
Note:	Chemical Formula: C ₁₄ H ₁₁ F ₃ N ₂ O ₂ Mol. Mass: 296.24 g/mol



[View online »](#)

Product images:

Flunixin



2-[[2-Methyl-3-(trifluoromethyl)phenyl]amino]pyridine-3-carboxylic acid