

## Product datasheet for **AP09733PU-L**

### Isoxsuprine Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.312µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Isoxsuprine (p-OH)-BTG
Specificity:	This antibody reacts to Isoxsuprine. <b>Cross-Reactivity;</b> Isoxsuprine 100% Ritodrine 0.1% Dobutamine 0.1% Fenoterol 0.1% Ractopamine 0.1%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Isoxsuprine is a synthetic chemical with epinephrine-Like effects. Its action relates specifically to the walls of arteries which are located in the muscles, internal organs, and extremities.
Synonyms:	Vasodilian, Duvadilan



[View online »](#)