

## Product datasheet for **AP09736PU-L**

### Kebuzone Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 0.625 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Kebuzone-BTG
Specificity:	Recognizes Kebuzone
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Kebuzone (or ketophenylbutazone) is a non-steroidal anti-inflammatory drug.
Synonyms:	ketophenylbutazone, Chetazolidin, Chebutan, Chepirol, Copirene, Hichillos, Ketazone, Recheton, Ketason, Chetil



[View online »](#)