

## Product datasheet for **AP09737PU-L**

### Ketoprofen Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 10 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Ketoprofen-BTG
Specificity:	This antibody recognizes Ketoprofen.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: One year from despatch.
Background:	Ketoprofen, (RS)2-(3-benzoylphenyl)-propionic acid, is a common non-steroidal anti-inflammatory drug (NSAID) with analgesic and antipyretic effects. It acts by inhibiting the body's production of prostaglandin.
Synonyms:	Ketoprofene, Profenid, Alrheumun, Capisten, Profenid, Oruvail, Actron, Epatec



[View online »](#)