

## Product datasheet for **AP09739PU-L**

### Lasalocid A Sheep Polyclonal Antibody

#### Product data:

|                       |   |
|-----------------------|---|
| Product Type:         | Primary Antibodies  |
| Applications:         | ELISA   |
| Recommended Dilution: | ELISA: 5µg/ml.  |
| Host:                 | Sheep   |
| Isotype:              | IgG   |
| Clonality:            | Polyclonal  |
| Immunogen:            | Lasalocid A-BTG   |
| Specificity:          | This antibody reacts to Lasalocid A.  |
| Formulation:          | 20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative<br>State: Ig Fraction<br>State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures |
| Concentration:        | lot specific  |
| Conjugation:          | Unconjugated  |
| Storage:              | Store the antibody at -20°C.<br>Avoid repeated freezing and thawing.  |
| Stability:            | Shelf life: one year from despatch.   |



[View online »](#)