

Product datasheet for **AP09746PU-L**

Megestrol Acetate Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 2.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Megestrol acetate (3)-BTG
Specificity:	This antibody reacts to Megestrol Acetate.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Megesterol: a synthetic progestational compound used to treat endometrial carcinoma and an appetite stimulant approved for the treatment of weight loss in people with AIDS. Megestrol acetate is a synthetic version of the female hormone PROGESTERONE. Most of the weight gain it leads to has been found to be fat rather than PROTEIN. Possible SIDE EFFECTS include HYPERGLYCEMIA, occasional impotence and decreased LIBIDO in men and uterine bleeding in women. Megestrol acetate is a synthetic progestogen used in the palliative management of recurrent, inoperable or metastatic endometrial carcinoma.



[View online »](#)