

Product datasheet for **AP09753PU-L**

Mescaline Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 2.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Mescaline-BTG
Specificity:	This antibody reacts to Mescaline.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Mescaline is a naturally occurring psychedelic alkaloid of the phenethylamine class. It is mainly obtained from members of the Cactaceae plant family and is also found in small amounts in certain members of the Fabaceae (bean) family. Users typically experience visual hallucinations and radically altered states of consciousness, often experienced as pleasurable and illuminating but occasionally is accompanied by feelings of anxiety or revulsion. Other effects include: open and closed eye visualizations, euphoria, dream-like state, laughter and a psychedelic experience.
Synonyms:	3,4,5-trimethoxyphenethylamine



[View online »](#)