

Product datasheet for **AP09756PU-L**

Methadone Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.08 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Methadone(p)-BSA
Specificity:	Recognizes Methadone
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Methadone is a rigorously well-tested medication that is safe and efficacious for the treatment of narcotic withdrawal and dependence. For more than 30 years this synthetic narcotic has been used to treat opioid addiction. Heroin releases an excess of dopamine in the body and causes users to need an opiate continuously occupying the opioid receptor in the brain. Methadone occupies this receptor and is the stabilizing factor that permits addicts on methadone to change their behavior and to discontinue heroin use.



[View online »](#)