

Product datasheet for **AP09762PU-L**

Methaqualone (4') Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.15µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Methaqualone (4')-BSA
Specificity:	This antibody reacts to Methaqualone (4').
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	Methaqualone is a nonbarbiturate sedative/hypnotic that produces sleep for about 6 to 8 hours. It also produces muscular relaxation, feelings of contentment, and total passivity. It is an addictive drug.



[View online »](#)