

## Product datasheet for **AP09787PU-N**

### Oxazepam Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 5.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Oxazepam (3)-BSA
Specificity:	This antibody reacts to Oxazepam (3). <b>Cross-Reactivity:</b> Oxazepam 100% Lorazepam 0.77% Temazepam 0.77% Alprazolam 0.77% alpha-Hydroxyalprazolam 0.77% Nordiazepam 974.68% Nitrazepam 55% 7-Aminonitrazepam 0.77% Clonazepam 0.77%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Benzodiazepine, BZO



[View online »](#)