

Product datasheet for **AP09796PU-L**

Phenobarbital (m/p) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 6.74µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Phenobarbital (m/p)-BSA
Specificity:	This antibody reacts to Phenobarbital (m/p). Cross-Reactivity: Phenobarbital 100% Secobarbital 2% Pentobarbital 8% Butalbital 5% Amobarbital 7% Barbital 5%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Luminal, Phenobarbitone, 5-ethyl-5-phenylpyrimidine-2,4,6(1H,3H,5H)-trione



[View online »](#)