

Product datasheet for **AP09800PU-L**

Progesterone (11-BSA) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.0 µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Progesterone (11)-BSA
Specificity:	This antibody recognizes Progesterone (11).
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	20-dione, pregn-4-ene-3



[View online »](#)