

## Product datasheet for **AP09804PU-N**

### Ractopamine (10') Sheep Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 0.625µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Ractopamine (10')-BSA
Specificity:	This antibody reacts to Ractopamine (10'). <b>Cross-Reactivity:</b> Ractopamine 100% Dobutamine 0.8% Isoxsuprine 0.6% DL-Isoproterenol 0.5% Ritodrine 0.5% Salmeterol Xinafoate 0.5% Methyl Clenbuterol 0.05% Fenoterol 0.05% Pirbuterol 0.05% Tulobuterol 0.05% Dopamine 0.01% Salbutamol 0.01% Clenbuterol 0.01% Metaproterenol 0.01% Terbutaline 0.01%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific



[View online »](#)

<b>Conjugation:</b>	Unconjugated
<b>Storage:</b>	Store the antibody at -20°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.