

Product datasheet for **AP09805PU-N**

Ractopamine (10) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 1.25µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Ractopamine (10)-BSA
Specificity:	This antibody reacts to Ractopamine (10). Cross-Reactivity: Ractopamine 100% Dobutamine 0.8% Isoxsuprine 0.6% DL-Isoproterenol 0.5% Ritodrine 0.5% Salmeterol Xinafoate 0.5% Methyl Clenbuterol 0.05% Fenoterol 0.05% Pirbuterol 0.05% Tulobuterol 0.05% Dopamine 0.01% Salbutamol 0.01% Clenbuterol 0.01% Metaproterenol 0.01% Terbutaline 0.01%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific



[View online »](#)

Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.