

Product datasheet for **AP09812PU-N**

Secobarbital Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 5.75µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Secobarbital (5)-BSA
Specificity:	This antibody reacts to Secobarbital. Cross-Reactivity: Secobarbital 100% Phenobarbital 23.1% Pentobarbital 23.1% Butalbital 23.1% Amobarbital 23.1% Barbital 23.1%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)