

Product datasheet for **AP09815PU-L**

Stanozolol Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 6.6µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Stanozolol (17)-BSA
Specificity:	This antibody reacts to Stanozolol.
Cross-Reactivity;	
	Stanozolol 100%
	16 beta-Hydroxystanozolol 40%
	4 alpha-Hydroxystanozolol 0.9%
	3-Hydroxy-17-Epi-Stanozolol 0.1%
	3-Hydroxystanozolol 0.1%
	19-Nortestosterone (17alpha) 0.1%
	19-Nortestosterone (17beta) 0.1%
	Trenbolone (17alpha) 0.1%
	Trenbolone (17beta) 0.1%
	19-Nor-4-Androstene - 3,17-Dione 0.1%
	Boldenone 0.1%
	Progesterone 0.1%
	Testosterone (17beta) 0.1%
	Cortisol 0.1%
	Estriol 0.1%
	Estrone 0.1%
	Danazol 0.1%
	Methandriol 0.1%
	17alpha-Methyltestosterone 0.1%



[View online »](#)

Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Winstrol, Androstanazole, Methylstanazole