

Product datasheet for **AP09829PU-L**

Sulphaphenazole Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 10.0 µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Sulfaphenazole-BTG
Specificity:	This antibody reacts to Sulphaphenazole.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Sulfaphenazole, Sulfabid, Sulfaphenazon, Raziosulfa, Firmazolo, Paidazolo, Plisulfan, Depocid



[View online »](#)