

Product datasheet for **AP09831PU-L**

Sulphaquinoxaline Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 10µg/ml
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Sulphaquinoxaline-HSA
Specificity:	This antibody reacts to Sulphaquinoxaline. Cross-Reactivity: Sulphaquinoxaline 100% Sulphamethazine 0.01% Sulphamerazine 0.01% Sulphadimethoxine 0.01% Sulphadiazine 0.01% Sulphapyridine 0.01% Sulphaguanidine 0.01% Sulphanilamide 0.01% Sulphamethizole 0.01% Sulphachloropyridazine 0.01% Sulphamethoxypyridazine 0.01% Sulphadoxine 0.01% Sulphisoxazole 0.01%
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated



[View online »](#)

Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Sulfaquinoxaline, Avicoid, Sulquin, Sulfabenzpyrazine, Ursokoxaline, Italquina, Kokozigal, Sulfaline, Sulfacox, Sulphaquinoxaline