

Product datasheet for **AP09836PU-L**

Temazepam Sheep Polyclonal Antibody

Product data:

| | |
|-----------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA |
| Recommended Dilution: | ELISA: 3.5µg/ml. |
| Host: | Sheep |
| Isotype: | IgG |
| Clonality: | Polyclonal |
| Immunogen: | Temazepam (3)-BSA |
| Specificity: | This antibody reacts to Temazepam. |
| Formulation: | 20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures |
| Concentration: | lot specific |
| Conjugation: | Unconjugated |
| Storage: | Store the antibody at -20°C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| Synonyms: | Restoril, Hydroxydiazepam, Methyloxazepam, Oxydiazepam, Crisonar, Euhypnos, Levanxol, Signopam, Levanxene, Mabertin, N-Methyloxazepam, Dasuen, Nortem , Planum, Temaz, Tenox |



[View online »](#)