

Product datasheet for **AP09836PU-N**

Temazepam Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	ELISA: 3.5µg/ml.
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	Temazepam (3)-BSA
Specificity:	This antibody reacts to Temazepam.
Formulation:	20mM Phosphate, 150mM Sodium Chloride, pH 7.2 containing 0.09% Sodium Azide as preservative State: Ig Fraction State: Liquid Ig fraction prepared by Caprylic Acid and Ammonium Sulphate precipitation procedures
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	Restoril, Hydroxydiazepam, Methyloxazepam, Oxydiazepam, Crisonar, Euhypnos, Levanxol, Signopam, Levanxene, Mabertin, N-Methyloxazepam, Dasuen, Nortem , Planum, Temaz, Tenox



[View online »](#)